### **Endless possibilities**

Take advantage of iOptima INT's extreme versatility and customise your system according to your needs. For more information, please contact Bien-Air USA on (800) 433-2436.

Kits soon available:







## DISCOVER THE BEST OF **SWISS TECHNOLOGY**













SWISS 🛨 MADE



# **INEW IOPTIMA**

WaveOne® READY









**ENDO DONTICS** WaveOne®



The iOptima INT electric system is incredibly versatile and enables you to perform all your dental operations - endodontics (including WaveOne®), restoration, oral surgery, implantology and scaling. Infinitely practical, the iOptima INT can be easily installed on any existing dental unit and can accommodate up to 2 motors and one Piezo scaler.

SWISS # MADE

www.bienair.com

limitless versatility

### All in One Make your life easier with Bien-Air Dental's all-in-one system. Highly versatile, the iOptima INT enables you to perform most of your dental operations. **ENDO RESTO IMPLANT SCALING** DONTICS PIEZ0 WaveOne 03.0 Fully integrated into your dental unit, the iOptima INT system guarantees optimal space savings and ensures that your working area is fully operational. The iOptima INT is an all-in-one space-saving

#### **Endodontics**

Coupled to the iOptima INT and MX2 motor, the new Bien-Air Dental CA ENDO handpiece has been specially designed for root canal treatment in continuous as well as WaveOne® mode. It is therefore perfectly compatible with standard endodontic files as well as with WaveOne® files by Dentsply Sirona. In addition, its small head ensures optimal handling and improved visibility of the working area.



100%

S ()

. . . . . . . . . .



Make your root canal treatment procedures easier with the WaveOne® mode by Dentsply Sirona. Thanks to alternating movements in the cutting direction and reverse direction, the file moves smoothly along the root canal and reduces the risk of file breakage. Your endodontic treatment are therefore safer and more successful.

The WaveOne® endodontics mode is available as an option.

NEW CA ENDO Bien-Air



and money-saving solution.

#### Restoration

Thanks to its multiple advantages, the iOptima INT, MX2 motor and EVO.15 1:5 highspeed handpiece combination allows you to carry out all your restoration work, including the treatment of caries and finishing / polishing work.

With advanced electronic management, the MX2 is extremely stable and powerful, even at low speed. Combined with the EVO.15 1:5's revolutionary file clamping system, this motor/ handpiece ensures vibration-free operation with unmatched accuracy.

# Bien Air iOptima WaveOne® (RECIPROC® ENDODONTICS 20'000 \*

Built directly into the dental unit, the iOptima INT system can be controlled either using an iPod touch or using an iPad mini. The refined ergonomics and user-friendliness of these tablets will enable you to focus on what matters most: the patient.

In addition, the numerous pre-programmed customisable sequences of the iOptima INT application will make your treatments considerably easier. Let us guide you!



#### Oral surgery & implantology

Driven by the iOptima<sup>INT</sup>, the MX-i motor and CA 20:1 handpiece combination offers you the best technology for all your implantology and oral surgery procedures.

While the MX-i's high torque will overcome the most critical implant operations, the CA 20:1's internal irrigation system will help you avoid discomfort when using the instrument.











